


## Aspartate Aminotransferase (AST)

### Test principle: UV test



Method according to the International Federation of Clinical Chemistry (IFCC), without pyridoxal-5'-phosphate. AST in the sample catalyzes the transfer of an amino group between L-aspartate and  $\alpha$ -ketoglutarate to form oxaloacetate and L-glutamate. The oxaloacetate then reacts with NADH in the presence of malate dehydrogenase (MDH) to form  $\text{NAD}^+$ .

The rate of the NADH oxidation is directly proportional to the catalytic AST activity. It is determined by measuring the decrease in absorbance at 340 nm.

### Proposed reagent composition

approximately 2+1 formulation

#### Reagent 1

Composition	Concentration	Catalog Number
Buffer (TRIS, pH 7.8)	280 mmol/l	10 153 265 001
L-Aspartate	800 mmol/l	
Malate dehydrogenase (MDH)	>2.1 kU/l	10 267 155 103
Lactate dehydrogenase (LDH)	>3.0 kU/l	10 679 666 103 or 12 235 650 103
Detergent, preservative, stabilizer, such as Sodium azide		
Triton X-100	0.01 %	10 743 119 103
Albumin	0.25 %	10 738 328 103

#### Reagent 2

Composition	Concentration	Catalog Number
NADH	>1.7 mmol/l	10 004 642 103
$\alpha$ -Ketoglutarate, di-Na	94 mmol/l	10 040 584 103
Preservative, such as Sodium azide		

Products are for further processing only.

All brands or product names are trademarks of their respective holders.

[custombiotech.roche.com](http://custombiotech.roche.com)

## Your Roche Custom Biotech Customer Service

### **Europe, Middle East, Africa, Latin America**

Phone +49 621 759 8580

Fax +49 621 759 8610

[mannheim.custombiotech@roche.com](mailto:mannheim.custombiotech@roche.com)

**Japan** Phone +81 3 5443 5285

Fax +81 3 5443 7934

[japan.custombiotech@roche.com](mailto:japan.custombiotech@roche.com)

**Asia Pacific** Phone +65 6371 6638

Fax +65 6371 6601

[apac.custombiotech@roche.com](mailto:apac.custombiotech@roche.com)

### **United States**

Phone +1 800 428 5433, ext. 14649 (toll-free)

Fax +1 317 521 4065

[custombiotech.ussales@roche.com](mailto:custombiotech.ussales@roche.com)

**Canada** Phone +1 450 686 7050

Fax +1 450 686 7012

[custombiotech.can@roche.com](mailto:custombiotech.can@roche.com)

### **Published by**

Roche Diagnostics GmbH

Sandhofer Straße 116

68305 Mannheim

Germany

© 2011 Roche Diagnostics.

All rights reserved.

05837723990 © 0311