


HawkZ05 Fast Polymerase

Single enzyme for fast and efficient RNA & DNA amplification

The choice of polymerase is a crucial component to success in amplification, for both RNA and DNA. CustomBiotech offers a single-enzyme solution for Real-Time PCR detection of both targets – HawkZ05 Fast Polymerase. HawkZ05 Fast Polymerase enables you to run fast one-step RT-PCR protocols with just one enzyme. The enzyme aptamer-oligonucleotide mixture is a reversible, temperature-

dependent hot start system. The aptamer acts as a molecular switch, changing its temperature-dependent tertiary structure. Dropping the temperature below +55°C shuts off the polymerase activity, while temperatures above +60°C fully activate the enzyme.



HawkZ05 Fast Polymerase benefits:

- **Achieve improved convenience and flexibility in IVD PCR workflows**

HawkZ05 Fast Polymerase hot start system enables simultaneous detection of both RNA and DNA targets, such as HCV and HPV.

- **Develop and run assays for difficult targets**

Achieve reliable amplification of your low-copy and complex RNA targets due to high temperature reverse transcription at 60–65°C and improved RNA processivity.

- **Reduce time to result**

Use a fast PCR protocol to reduce your total assay and turnaround time (TAT) without losing sensitivity and linearity. Run RT-PCR reaction in as little as 44 minutes.

- **Easily work with crude samples**

Hawk Z05 Fast Polymerase shows high inhibitor tolerance towards PCR inhibitors such as Heparin and SDS, and robust performance in problematic clinical samples from, for example, stool or nasal swabs.

- **Suitable for various diagnostic PCR systems and applications**

Glycerol-free enzyme available for lyophilization and improved liquid handling.

- **Obtain consistent results**

Rely on our standardized manufacturing processes, including extensive QC release testing and thus, high lot-to-lot consistency, providing the perfect basis for (IVD) kit manufacturers and end users.

Ordering information

Product	Pack size	Catalog number
HawkZ05 Fast DNA Polymerase, 40 U/µl	Custom fill	07731264103
HawkZ05 Fast DNA Polymerase, Glycerol-free, 200 U/µl	Custom fill	07731329103

Related products	Pack size	Catalog number
NucleoMixes PCR Grade		
dATP, dCTP, dGTP, dTTP (25 mmol/l each)	20 ml 100 ml	04920171103 03991016103
dATP, dCTP, dGTP, dTTP (10 mmol/l each)	100 ml	03186083103
dATP, dCTP, dGTP, dUTP (25 mmol/l each)	20 ml	04980905103
dATP, dCTP, dGTP, dUTP (10 mmol/l each)	100 ml	03186075103

Regulatory disclaimer

For further processing only.

Trademark

HAWKZ05 is a trademark of Roche.

custombiotech.roche.com

Please contact your local CustomBiotech representative

Europe, Middle East, Africa, Latin America

Phone +49 621 759 8580
Fax +49 621 759 6385
mannheim.custombiotech@roche.com

United States

Phone +1 800 428 5433, ext. 14649 (toll-free)
Fax +1 317 521 4065
custombiotech.ussales@roche.com

Canada

Phone +1 450 686 7050
Fax +1 450 686 7012
custombiotech.can@roche.com

Japan

Phone +81 3 5443 5285
Fax +81 3 5443 7934
japan.custombiotech@roche.com

Asia Pacific

Phone +65 6371 6638
Fax +65 6371 6601
apac.custombiotech@roche.com

Published by

Roche Diagnostics GmbH
Sandhofer Straße 116
68305 Mannheim
Germany

© 2016 Roche Diagnostics GmbH
All rights reserved.